




APHTGEL
RemeSense

PRODUCT:

Aphtgel RemeSense for mouth ulcers, is the new gel solution for immediate pain relief and treatment of mouth ulcers, cheek bites, abrasion from braces and denture irritation.

The gel has active properties that help to heal (recurrent) mouth ulcers and oral irritations.

Apply the gel to the affected area and it will form a protective barrier over the mouth ulcer while the triple action formula helps to speed up the natural healing. The gel is very easy to use, which makes it extremely suitable for children and elderly people. Aphtgel RemeSense works fast and is effective.

TECHNOLOGY:

Aphtgel RemeSense is based on Sylphar's film forming technology. Upon application, it will form a transparent, thin muco-adhesive film on the mouth ulcer. This film will protect and isolate the injured mouth area affected by the ulcer.

The main components of the device are hyaluronate macromolecules and film forming polymers. Hyaluronate macromolecules possess an extremely high capacity to bind water and are able to form a viscous, protective macro-aggregate over the affected tissues. Based on those properties, hyaluronate macromolecules accelerate the rate of tissue reconstruction of damaged mouth area.

Other important components are film forming polymers: Calcium/Sodium PVM/MA Copolymer, Calcium Carboxymethyl Cellulose and Xanthan gum:

- PVM/MA Copolymers are widely used in pharmaceutical applications for their excellent film forming properties and high quality bio-adhesive performances.
- Cellulose is an excellent film former and in the field of wound and tissue regeneration cellulose is being used to stimulate wound healing.
- Xanthan gum is hetero-polysaccharide with a very high molecular weight. That's why they are unable to penetrate the skin and will form a protective barrier on the skin..

FEATURES:

- ✓ Easy to use
- ✓ Reduces pain of mouth ulcers and oral irritations
- ✓ Stimulates healing
- ✓ Rapid, soothing effect

CONTENT:

1 x 10ml Aphtgel

Sylphar

Sylphar N.V.
Xavier De Cocklaan 42
B-9831 Deurle, Belgium
Tel +32 9 321 70 80
Fax +32 9 321 70 90
info@sylphar.com
www.sylphar.com
www.afterremesense.com